

PRO-MULTI USE BOOT ORTHOSIS Fitting INSTRUCTIONS

Fitting Instructions: LG PRO-MULTI USE BOOT ORTHOSIS. Adult and Pediatric

Custom fit without tools or heat. Bend it as you see fit. Open all Velcro closures to Adjust. See Video on-line www.LEEDerGroup.com

Open all Velcro Closures and Adjust to Fit Condition. Hand Bend as needed without tools or heat.



Adult/Pediatric Sizes

LMUAAdult

LMUPediatric

Congratulations: you've chosen a product with Flex-E-Core technology: the Multi Use Boot. There is no other Adjustable boot product like this on the market. Custom fit without tools or heat. Bend to fit.

Product Codes: LMUAAdult, **LMUP**ediatric

Make the PRO-MULTI USE BOOT ORTHOSIS fit perfectly, BEND IT TO FIT. Custom fitting is EASY:

PRO-MULTI USE BOOT ORTHOSIS Fitting INSTRUCTIONS

1. If the FOOT section is too long or too short, disconnect the side (dorsiflexion) adjustment straps and adjust the point of bend above or below the heel to shorten or lengthen the splint.

2. Then, RE-BEND the boot so that the FOOT section fits the user, CUSTOM FIT it either longer or shorter.

3. Inversion or Eversion: twist the whole frame to match and correct the problem.

4. Night Splinting: wear with sock; adjust the dorsiflexion straps to counteract plantar flexion contracture; apply pressure in small increments. Cut off straps as desired. **OPTION:** use foam dorsiflexion wedge for extra toe dorsiflex.

5. For HEEL FLOAT make small bends along the “heel” surface of the boot. Make small bends and small changes to get the perfect fit.

6. For Best Bending Results: USE a “rounded” edge, such as the side rails of a bed or top of the backrest to a chair.

Cut off straps as needed.

**Mold Boot by
HAND as
Necessary to
Custom Fit:
Heel Float,
Plantar Fasciitis,
Night Splint,
Inversion,
Eversion**

New

**Softest Adjustable
Interface Available
Flex-E-Core
Technology**

Contraindications: at least EACH SHIFT check for swelling or redness; check each hour during first shift.

Contact your supplier to re-order.