

SEWHO INSTALLATION INSTRUCTIONS

"Since 1996"




Give More Care for LESS COST ... The LEEDer way Installing the Easy Fit SEWHO No Tools Needed * Flex-E-Core



1. Loosen the **YELLOW wing bolts** on the clamping device to slide the bottom plate under the wheelchair ARM REST. Remove one nut if necessary to get under.
2. Once the bottom plate passes under the wheelchair armrest, push the plate into the

Insert "teeth" of bottom plate into slots



protective, black foam to engage the two parts. Hook the plate as high as possible on the row of slots to secure the total assembly. Tighten the **YELLOW wing bolts** to secure the unit to the wheelchair. **IF NECESSARY**, use wood spacer block under wheelchair arm for better fit.

Connect bottom plate to side slots

3. With the patient in the wheelchair place his/her arm in the arm tray. Based upon the patient's needs, locate the proper position for the SEWHO on the wheelchair **ARM REST** by sliding it forward or backward. After selecting position, push downward on the clamping device while tightening the **YELLOW wing bolts**.

Short Arm Tray




4. **Flex-E-Core Arm Tray:** Firmly grasp the Arm Tray with one hand and the "Elbow Stopper" end of the arm tray with the other; rotate the "Elbow Stopper" toward the outside of the chair away from the patient to prevent the elbow from sliding out of the Tray. The Universal arm tray can be adjusted to fit either Right or Left.

NEW!
Allen Screw to control rotation

Adjust Arm Tray as necessary.

5. For Vertical Height adjustment remove the patient's arm from the arm tray. Loosen both **BLACK wing bolts** on the flat side bar 3 full turns.



6. Hold the Arm tray at the desired level and tighten the two **Black wing bolts**. "View Holes" are provided next to the wing bolts to help align the wing bolts with the holes. Tighten the wing bolts.

Vertical Elevation Adjustment



7. Adjust the **RED wing bolts** to control lateral and Pivot movement after the desired vertical position is reached. **TO LOCK** the Lateral-Pivot movement, loosen the **RED wing bolts** and tighten the small **Silver screws** with finger or Allen wrench provided; then, Re-tighten **RED wing bolts**.

Horizontal Adjustment allows rotation for THERAPY

ASSEMBLY WILL NOT support a patient's FULL body weight!

Toll Free Customer Service 866-814-0192
www.LEEDerGroup.com

Toll Free FAX 866-818-0373
LEEDerGroup@LEEDerGroup.com

8508 NW 66th St Miami FL 33166