

Custom Fitting INSTRUCTIONS

"Since 1996"



Give More Care for LESS COST ... The LEEDer way

Fitting

NEW & IMPROVED

MULTI USE
BOOT /
NIGHT SPLINT



New

Softest Adjustable
Interface Available
Flex-E-Core
Technology

Dorsiflexion
Adjustment
Straps. One
on each side.

Anti
Rotation
Bar.
Rotate for
Position

Foot Cover
Adjustment
Strap

Calf Strap
Adjustment

Walking Sole
Simply wipe clean.

Instructions: Multi Use Boot. LMUA, P
Custom fit without tools or heat. Bend it as you see fit.
Open all Velcro closures to Adjust.

Open all Velcro Closures and Adjust to Fit Condition. Hand Bend as needed without tools or heat.

Adult/Pediatric Sizes

LMUA^{adult}

LMUP^{pediatric}

Fitting Instructions: Multi Use Boot. LMUA, P

***Congratulations: you've chosen a product with
LEEDer's Flex-E-Core technology: the Multi Use Boot.***

Toll Free Customer Service 866-814-0192
www.leedergroup.com

Toll Free FAX 866-818-0373
leedergroup@leedergroup.com

8508 NW 66th St Miami FL 33166

Custom Fitting INSTRUCTIONS “Since 1996”

Give More Care for LESS COST ... The LEEDer way

There is no other Adjustable boot product like this on the market. Custom fit without tools or heat. Bend to fit.

Product Codes: **LMUA**adult, **LMUP**ediatric



Make the Multi Use Boot fit perfectly, BEND IT TO FIT. Custom fitting is EASY:

1. If the FOOT section is too long or too short, disconnect the side (dorsiflexion) adjustment straps and adjust the point of bend above or below the heel to shorten or lengthen the splint.

2. Then, RE-BEND the boot so that the FOOT section fits the user, CUSTOM FIT it either longer or shorter.

3. Inversion or Eversion: twist the whole frame to match and correct the problem.

4. **Night Splinting:** wear with sock; adjust the dorsiflexion straps to counteract plantar flexion contracture; apply pressure in small increments. Cut off straps as desired. **OPTION:** use foam dorsiflexion wedge for extra toe dorsiflex.

5. For HEEL FLOAT make small bends along the “heel” surface of the boot. Make small bends and small changes to get the perfect fit.

6. For Best Bending Results: USE a “rounded” edge, such as the side rails of a bed or top of the backrest to a chair.

Cut off straps as needed.