

CAM WALKER FITTING INSTRUCTIONS

Indications/Contradictions

Before applying the brace, please read these instructions completely and carefully. Correct application is vital to the proper functioning of the brace.

INDICATIONS

- Acute ankle sprains
 - Soft tissue injuries of the lower leg
 - Stress fractures of the lower leg
 - Stable fractures of the foot and ankle
 - Soft tissue injuries of the lower leg
 - Stable delayed union or non-union fractures of the distal tibia and fibula
- Achilles tendon repairs

CONTRADICTIONS

Application of this device is recommended only when the fracture is demonstrably stable and there are acceptable limits of angular and rotational deformity. The determination of when to apply the brace rests strictly at discretion of the treating physician.

This product is designed to complement the variety of medical treatments common to the above afflictions. The frequency and duration of use should be determined by your prescribing healthcare professional.



TOLL FREE: 866-814-0192

A

For first time application, loosen straps and remove liner from boots

B

Place foot in liner and secure with contact closure. Make sure the heel fits snugly into the posterior portion of the liner, again making sure that it is snug from bottom to top. Some trimming of the liner may be necessary to accomplish optimum fit.

C

Spread the uprights using both hands and step into boot, aligning uprights with midline of the ankle. **Note:** if uprights have plastic tube liner remove them now to affix upright to liner.

D

Secure boot straps starting at the toes and work up the leg.